

Laboquick

Urinalysis Reagent Strips

10 Parameter

Glucose

Bilirubin

Ketone

Specific Gravity

Blood

pH

Protein

Urobilinogen

Nitrite

Leukocytes

HAM.LUET.01

100 Strips

- For in vitro diagnostic use only.
- Read instructions for use before starting.
- Do not use after the expiration date.
- Keep device out of the reach of the children.
- Dispose used materials as biological waste.
- Keep between +2...+30°C. Do not freeze.
- Keep away from radiation, moisture, odour.
- Only take required amount of strips that will be used.
- The strip should remain in the closed canister until use.
- Do not touch the reagent areas of the strip.



KÖROĞLU

KÖROĞLU TIBBİ MALZ. ve KOZMETİK SAN. TİC. LTD. ŞTİ.
 1775/23 Sok. No:4/A Mevlana Mah. Bornova / İZMİR
 Tel: (0232) 388 73 74 • Faks: (0232) 388 40 43
 www.koroglugroup.com • info@koroglugroup.com

